


# ZEN-PX2

## Portable Dental X-ray System

- Excellent Image Quality
- Elegant Design
- Ideal Output for Film & Digital sensor : 60kV, 2mA
- Easy to use

### Compact Outside, Smart Inside

New era of Intra-Oral X-ray system is now in your hand.  
ZEN-PX2 will offer you ease and comfort in your practice.




New era of Intra-Oral X-ray system is now in your hand. ZEN-PX2 will offer you ease and comfort in your practice.

- Light and compact unit lets you bring your practice to anywhere.
- High performance Tube and HFG makes sharp and clear images.
- Advanced Radiation Shielding ensures less than 0.13% of allowed doses leak.
- Ideal to use with Digital sensor and chemical Film.

## Technical Specification & Dimension

<b>X-ray Generator</b>	Type	High Frequency Inverter
	Power Output	180W
<b>X-ray Installation</b>	Tube Voltage	60 kV
	Tube Current	2mA
	Heat Capacity	8.5 KHU
<b>X-ray Tube</b>	Anode Type	Stationary
	Anode Angle	20°
	Focal Spot	0.8mm
<b>X-ray Control</b>	Exposure Time Set	56 Steps (0.01 ~ 2.0 sec)
<b>Power Requirement</b>		22.2 VDC
<b>Inherent Filtration</b>	Al equivalent	1.8 mmAl
<b>Battery Capacity</b>	Number of exam	100 per full charge
<b>Weight</b>	Main Body	5.18lb (2.35kg)
<b>Dimension</b>	Main Body	7.75" x 5.78" x 5.70" (19.7cm x 14.7cm x 14.5cm)

\* Please do not hesitate to contact us if you need detailed specification of our products.



### Components of ZEN-PX2

#### Main Unit

- X-ray tube assembly,  
Beam limiting device,  
Battery, Grip tie

#### Accessories

- Power adapter, Neck strap,  
Remote control, Manual

#### ● Compact Design

Camera type ZEN-PX2 offers convenient diagnosis with its portable ability and compact design. ZEN-PX2 also offers large battery capacity capable of over 100 exams on single charge.

#### ● Graphic LCD

Large LCD display lets you monitor the process and setting.

#### ● Grip Tie

Adjust length of this comfortable tie to fit your hand & be completely mobile.

#### ● Operation Panel

X-ray operation have never been easier with easy to understand push button & optimized pre-settings.